

ROTOFIX 32A

BENCHTOP CENTRIFUGE

FOR **CLINICAL** & **RESEARCH** APPLICATIONS

FEATURES:

- Rugged, versatile, indispensable benchtop
- Multi-functional, general purpose centrifuge
- All-steel construction
- Run-Silent™ Operation (≤ 57 dB, softer than normal conversation)
- Easy-Lift™ / Safety-Lock™ Lid
- Five (5) Year Warranty
- FDA registered centrifuge manufacturer

Hettich
CENTRIFUGES

NEW

CENTRIFUGE PACKAGES

- BLOOD
- CELL CULTURE
- URINALYSIS
- CYTOLOGY

(See reverse for more info)



Find out more by visiting:

www.hettweb.com/rotofix-32a



MAX. CAPACITY:
4 x 100 mL

DIMENSIONS (HxWxD) 10.11" x 14.4" x 16.9"

TUBE TYPES MICROLITER, ROUND BOTTOM, BLOOD, URINE, CONICAL, SKIRTED, SCREW TOP, CHROME BATH, CYTOLOGY

ROTOFIX 32A CENTRIFUGE PACKAGES

FOR CLINICAL & RESEARCH APPLICATIONS

(BELOW IS A LIST OF THE **MOST POPULAR** CENTRIFUGE PACKAGES AVAILABLE FOR THIS UNIT. **ADDITIONAL CONFIGURATIONS** ARE AVAILABLE DEPENDING ON YOUR APPLICATION. PLEASE INQUIRE FOR MORE DETAILS)

ROTOFIX 32A BLOOD TUBE PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
	20	1.6 - 7	13 x 75/ 100	4,000	2,558	(1) Centrifuge
	16	4 - 10	16 x 75/ 100		2,522	(1) 90° Swing-out Rotor (4-place)
	Pkg No.:			32ABLOOD		
	Add Bio-containment Lids			(-BC)		
						(4) Buckets
						(4) Adapters for 13 x 75/100 mm
						(4) Adapters for 16 x 75/100 mm

ROTOFIX 32A CELL CULTURE PACKAGE - 4

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
	16	15 Conical	17 x 120	4,000	2,612	(1) Centrifuge
	4	50 Conical	29 x 115		2,576	(1) 90° Swing-out Rotor (4-place)
	Pkg No.:			32ACELLCULTURE4		
	Add Bio-containment Lids			(-BC)		
						(4) Buckets
						(4) Adapters for 15 mL
						(4) Adapters for 50 mL

ROTOFIX 32A CELL CULTURE PACKAGE - 6

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
	6	15 Conical	17 x 120	4,000	2,701	(1) Centrifuge
		50 Conical	29 x 115			(1) 90° Swing-out Rotor (6-place)
	Pkg No.:			32ACELLCULTURE6		(6) Buckets
						(6) Adapters for 15 mL

ROTOFIX 32A URINALYSIS PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	G-FORCE	PACKAGE CONTENTS
	12	KOVA (12)	17 x 102	4,000	2,451	1) Centrifuge (1) 90° Swing-out Rotor (4-place) (4) Buckets (4) Adapters
	Pkg No.:		32AURINALYSIS			



Find out more by visiting:
www.hettweb.com/rotofix-32a

DISTRIBUTED BY:

Hettich
CENTRIFUGES

HETTICH NORTH AMERICA
Toll Free: 1-866-370-4388
Email: info@hettweb.com

→ Discover more solutions online at: www.hettweb.com/centrifuges

Copyright © 2012 Hettich Instruments, LP. All rights reserved. For the most up to date information on Hettich please visit hettweb.com.
Printed in the USA using soy based inks on FSC certified recycled paper (25% pre-consumer waste & 25% post consumer waste).

PS-ROTOFIX-32A-E RevA 01-2013